





DCX Antibody

Product Code	CSB-PA007687
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43602
Immunogen	Synthesized peptide derived from Human Doublecortin around the non-phosphorylation site of S378.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DCX; DBCN; LISX; Neuronal migration protein doublecortin; Doublin; Lissencephalin-X; Lis-X
Immunogen Species	Homo sapiens (Human)
Target Names	DCX

Image



Western Blot analysis of RAW264.7 cells using Doublecortin Polyclonal Antibody



Western Blot analysis of 293T using Doublecortin Polyclonal Antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.