





DCX Antibody

Product Code	CSB-PA007687
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43602
Immunogen	Synthesized peptide derived from Human Doublecortin around the non-phosphorylation site of S378.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DCX; DBCN; LISX; Neuronal migration protein doublecortin; Doublin; Lissencephalin-X; Lis-X
Immunogen Species	Homo sapiens (Human)
Target Names	DCX

Image



Western Blot analysis of RAW264.7 cells using Doublecortin Polyclonal Antibody



Western Blot analysis of 293T using
Doublecortin Polyclonal Antibody.