




GPR113 Antibody

Product Code	CSB-PA007703
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZF5
Immunogen	Synthesized peptide derived from the N-terminal region of Human GPR113.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPR113; PGR23; Probable G-protein coupled receptor 113; G-protein coupled receptor PGR23
Immunogen Species	Homo sapiens (Human)
Target Names	ADGRF3

Image



Western Blot analysis of HeLa cells using GPR113 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.