




OLFML2A Antibody

Product Code	CSB-PA007764
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q68BL7
Immunogen	Synthesized peptide derived from the Internal region of Human Photomedin-1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OLFML2A; Olfactomedin-like protein 2A; Photomedin-1
Immunogen Species	Homo sapiens (Human)
Target Names	OLFML2A

Image



Western Blot analysis of RAW264.7 cells using Photomedin-1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.