




OLFML2A Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA007764 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q68BL7 |
| Immunogen | Synthesized peptide derived from the Internal region of Human Photomedin-1. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | OLFML2A; Olfactomedin-like protein 2A; Photomedin-1 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | OLFML2A |

Image



Western Blot analysis of RAW264.7 cells using Photomedin-1 Polyclonal Antibody