




DLG4 Antibody

Product Code	CSB-PA007805
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P78352
Immunogen	Synthesized peptide derived from Human PSD-95 around the non-phosphorylation site of S295.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DLG4; PSD95; Disks large homolog 4; Postsynaptic density protein 95; PSD-95; Synapse-associated protein 90; SAP-90; SAP90
Immunogen Species	Homo sapiens (Human)
Target Names	DLG4

Image



Western Blot analysis of 293T cells using PSD-95 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.