



DOK1 Antibody

Product Code	CSB-PA007817
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99704
Immunogen	Synthesized peptide derived from Human Dok-1 around the non-phosphorylation site of Y398.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DOK1; Docking protein 1; Downstream of tyrosine kinase 1; p62(dok); pp62
Immunogen Species	Homo sapiens (Human)
Target Names	DOK1