




APLNR Antibody

Product Code	CSB-PA007868
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35414
Immunogen	Synthesized peptide derived from the Internal region of Human APLNR.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	APLNR; AGTRL1; APJ; Apelin receptor; Angiotensin receptor-like 1; G-protein coupled receptor APJ; G-protein coupled receptor HG11
Immunogen Species	Homo sapiens (Human)
Target Names	APLNR

Image



Western Blot analysis of 293 cells using APLNR Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.