




Phospho-AHR (S36) Antibody

Product Code	CSB-PA007980
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35869
Immunogen	Synthesized peptide derived from Human Ah Receptor around the phosphorylation site of S36.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)

Image



Western Blot analysis of HepG2 cells using Phospho-Ah Receptor (S36) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.