



A2M Antibody

Product Code	CSB-PA008007
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01023
Immunogen	Synthesized peptide derived from the C-terminal region of Human Macroglobulin α -2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	A2M; CPAMD5; FWP007; Alpha-2-macroglobulin; Alpha-2-M; C3 and PZP-like alpha-2-macroglobulin domain-containing protein 5
Immunogen Species	Homo sapiens (Human)
Target Names	A2M
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.