




FOXD4L1 Antibody

Product Code	CSB-PA008024
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NU39
Immunogen	Synthesized peptide derived from the C-terminal region of Human FoxD4L1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FOXD4L1; Forkhead box protein D4-like 1; FOXD4-like 1
Immunogen Species	Homo sapiens (Human)
Target Names	FOXD4L1

Image



Western Blot analysis of HuvEc cells using FoxD4L1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.