



MARK2 Antibody

Product Code	CSB-PA008036
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7KZI7
Immunogen	Synthesized peptide derived from Human MARK2 around the non-phosphorylation site of T596.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, IF, ELISA; IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MARK2; EMK1; Serine/threonine-protein kinase MARK2; ELKL motif kinase 1; EMK-1; MAP/microtubule affinity-regulating kinase 2; PAR1 homolog; PAR1 homolog b; Par-1b; Par1b
Immunogen Species	Homo sapiens (Human)
Target Names	MARK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.