




CTTN Antibody

Product Code	CSB-PA008046
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14247
Immunogen	Synthesized peptide derived from Human Cortactin around the non-phosphorylation site of Y446.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CTTN; EMS1; Src substrate cortactin; Amplaxin; Oncogene EMS1
Immunogen Species	Homo sapiens (Human)
Target Names	CTTN

Image



Western Blot analysis of HeLa cells using Cortactin Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.