



EPB41 Antibody

Product Code	CSB-PA008055
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11171
Immunogen	Synthesized peptide derived from Human 4.1R around the non-phosphorylation site of Y660.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EPB41; E41P; Protein 4.1; P4.1; 4.1R; Band 4.1; EPB4.1
Immunogen Species	Homo sapiens (Human)
Target Names	EPB41
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.