



Phospho-DMTN (S403) Antibody

Product Code	CSB-PA008056
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08495
Immunogen	Synthesized peptide derived from Human Dematin around the phosphorylation site of S403.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EPB49; DMT; Dematin; Erythrocyte membrane protein band 4.9
Immunogen Species	Homo sapiens (Human)
Target Names	DMTN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.