




Phospho-EPHA4 (Y596) Antibody

Product Code	CSB-PA008062
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54764
Immunogen	Synthesized peptide derived from Human EphA4 around the phosphorylation site of Y596.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EPHA4; HEK8; SEK; TYRO1; Ephrin type-A receptor 4; EPH-like kinase 8; EK8; hEK8; Tyrosine-protein kinase TYRO1; Tyrosine-protein kinase receptor SEK
Immunogen Species	Homo sapiens (Human)
Target Names	EPHA4

Image



Western Blot analysis of JK cells using Phospho-EphA4 (Y596) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.