




EPHB4 Antibody

Product Code	CSB-PA008071
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54760
Immunogen	Synthesized peptide derived from the Internal region of Human EphB4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EPHB4; HTK; MYK1; TYRO11; Ephrin type-B receptor 4; Hepatoma transmembrane kinase; Tyrosine-protein kinase TYRO11
Immunogen Species	Homo sapiens (Human)
Target Names	EPHB4

Image



Western Blot analysis of HY926 cells using EphB4 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.