




EPOR Antibody

Product Code	CSB-PA008074
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19235
Immunogen	Synthesized peptide derived from Human EpoR around the non-phosphorylation site of Y426.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EPOR; Erythropoietin receptor; EPO-R
Immunogen Species	Homo sapiens (Human)
Target Names	EPOR

Image



Western Blot analysis of JK cells using EpoR Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.