




Phospho-EPS15 (Y849) Antibody

Product Code	CSB-PA008076
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42566
Immunogen	Synthesized peptide derived from Human Eps15 around the phosphorylation site of Y849.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EPS15; AF1P; Epidermal growth factor receptor substrate 15; Protein Eps15; Protein AF-1p
Immunogen Species	Homo sapiens (Human)
Target Names	EPS15

Image



Western Blot analysis of JK cells using Phospho-Eps15 (Y849) Polyclonal Antibody