



ERBB2 Antibody

Product Code	CSB-PA008089
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04626
Immunogen	Synthesized peptide derived from Human Neu around the non-phosphorylation site of Y1139.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ERBB2; HER2; MLN19; NEU; NGL; Receptor tyrosine-protein kinase erbB-2; Metastatic lymph node gene 19 protein; MLN 19; Proto-oncogene Neu; Proto-oncogene c-ErbB-2; Tyrosine kinase-type cell surface receptor HER2; p185erbB2; CD antigen CD340
Immunogen Species	Homo sapiens (Human)
Target Names	ERBB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.