



Phospho-AKT1/AKT2/AKT3 (Y315) Antibody

Product Code	CSB-PA008130
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31749/P31751/Q9Y243
Immunogen	Synthesized peptide derived from Human Akt around the phosphorylation site of Y315.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AKT1; PKB; RAC; RAC-alpha serine/threonine-protein kinase; Protein kinase B; PKB; Protein kinase B alpha; PKB alpha; Proto-oncogene c-Akt; RAC-PK-alpha; AKT2; RAC-beta serine/threonine-protein kinase; Protein kinase Akt-2; Protein kinase B
Immunogen Species	Homo sapiens (Human)
Target Names	AKT1/AKT2/AKT3