



# Phospho-ERF (T526) Antibody

Product Code	CSB-PA008135
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50548
Immunogen	Synthesized peptide derived from Human ERF around the phosphorylation site of T526.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ERF; ETS domain-containing transcription factor ERF; Ets2 repressor factor; PE-2
Immunogen Species	Homo sapiens (Human)
Target Names	ERF