



Phospho-ESR1 (S305) Antibody

Product Code	CSB-PA008146
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03372
Immunogen	Synthesized peptide derived from Human ER α around the phosphorylation site of S305.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ESR1; ESR; NR3A1; Estrogen receptor; ER; ER-alpha; Estradiol receptor; Nuclear receptor subfamily 3 group A member 1
Immunogen Species	Homo sapiens (Human)
Target Names	ESR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.