



Phospho-ESR1 (S305) Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA008146 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P03372 |
| Immunogen | Synthesized peptide derived from Human ER α around the phosphorylation site of S305. |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | IF,ELISA;IF:1/200-1/1000.ELISA:1/20000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | ESR1; ESR; NR3A1; Estrogen receptor; ER; ER-alpha; Estradiol receptor; Nuclear receptor subfamily 3 group A member 1 |
| Species | Homo sapiens (Human) |
| Target Names | ESR1 |