




ETV6 Antibody

Product Code	CSB-PA008157
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41212
Immunogen	Synthesized peptide derived from Human TEL around the non-phosphorylation site of S257.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ETV6; TEL; TEL1; Transcription factor ETV6; ETS translocation variant 6; ETS-related protein Tel1; Tel
Immunogen Species	Homo sapiens (Human)
Target Names	ETV6

Image



Western Blot analysis of HuvEc cells using TEL Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.