




Cleaved-F2R (S42) Antibody

Product Code	CSB-PA008161
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25116
Immunogen	Synthesized peptide derived from the N-terminal region of Human Thrombin R.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	F2R; CF2R; PAR1; TR; Proteinase-activated receptor 1; PAR-1; Coagulation factor II receptor; Thrombin receptor
Immunogen Species	Homo sapiens (Human)
Target Names	F2R

Image



Western Blot analysis of Jurkat cells using Cleaved-Thrombin R (S42) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.