



Phospho-FANCD2 (S222) Antibody

Product Code	CSB-PA008177
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXW9
Immunogen	Synthesized peptide derived from Human FANCD2 around the phosphorylation site of S222.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FANCD2; FACD; Fanconi anemia group D2 protein; Protein FACD2
Immunogen Species	Homo sapiens (Human)
Target Names	FANCD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.