



OR8U1/OR8U8/OR8U9 Antibody

Product Code	CSB-PA008193
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NH10/P0C7N1/P0C7N5
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 8U1/8/9.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR8U1; Olfactory receptor 8U1; OR8U8; Olfactory receptor 8U8; OR8U9; Olfactory receptor 8U9
Immunogen Species	Homo sapiens (Human)
Target Names	OR8U1/OR8U8/OR8U9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.