



OR5L1/OR5L2 Antibody

Product Code	CSB-PA008195
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGL2/Q8NGL0
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 5L1/2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	IF, ELISA; IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR5L1; Olfactory receptor 5L1; OST262; Olfactory receptor OR11-151; OR5L2; Olfactory receptor 5L2; HTPCRX16; Olfactory receptor OR11-153
Immunogen Species	Homo sapiens (Human)
Target Names	OR5L1/OR5L2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.