




OR7E5P Antibody

Product Code	CSB-PA008196
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96N54
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 7E5P.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	Seven transmembrane helix receptor
Immunogen Species	Homo sapiens (Human)
Target Names	OR7E5P

Image



Western Blot analysis of 3T3 cells using Olfactory receptor 7E5P Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.