




OR8K3 Antibody

Product Code	CSB-PA008198
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NH51
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 8K3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR8K3; Olfactory receptor 8K3; Olfactory receptor OR11-181
Immunogen Species	Homo sapiens (Human)
Target Names	OR8K3

Image



Western Blot analysis of LOVO cells using Olfactory receptor 8K3 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.