




OR8J1 Antibody

Product Code	CSB-PA008199
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGP2
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 8J1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR8J1; Olfactory receptor 8J1; Olfactory receptor OR11-183
Immunogen Species	Homo sapiens (Human)
Target Names	OR8J1

Image



Western Blot analysis of Jurkat cells using Olfactory receptor 8J1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.