




OR4D6 Antibody

Product Code	CSB-PA008210
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGJ1
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 4D6.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR4D6; Olfactory receptor 4D6; Olfactory receptor OR11-250
Immunogen Species	Homo sapiens (Human)
Target Names	OR4D6

Image



Western Blot analysis of LOVO cells using Olfactory receptor 4D6 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.