




GPR115 Antibody

Product Code	CSB-PA008221
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZF3
Immunogen	Synthesized peptide derived from the C-terminal region of Human GPR115.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPR115; PGR18; Probable G-protein coupled receptor 115; G-protein coupled receptor PGR18
Immunogen Species	Homo sapiens (Human)
Target Names	ADGRF4

Image



Western Blot analysis of RAW264.7 cells using GPR115 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.