




ZSCAN25 Antibody

Product Code	CSB-PA008224
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NSZ9
Immunogen	Synthesized peptide derived from the Internal region of Human ZNF498.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZNF498; ZSCAN25; Zinc finger protein 498; Zinc finger and SCAN domain-containing protein 25
Immunogen Species	Homo sapiens (Human)
Target Names	ZSCAN25

Image



Western Blot analysis of A549 cells using ZNF498 Polyclonal Antibody