




CCZ1/CCZ1B Antibody

Product Code	CSB-PA008230
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P86791/P86790
Immunogen	Synthesized peptide derived from the C-terminal region of Human CGI-43.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CCZ1; C7orf28A; CGI-43; Vacuolar fusion protein CCZ1 homolog; CCZ1B; C7orf28B; Vacuolar fusion protein CCZ1 homolog B; Vacuolar fusion protein CCZ1 homolog-like
Immunogen Species	Homo sapiens (Human)
Target Names	CCZ1/CCZ1B

Image



Western Blot analysis of Jurkat cells using CGI-43 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.