




GPRC6A Antibody

Product Code	CSB-PA008232
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T6X5
Immunogen	Synthesized peptide derived from the Internal region of Human GPRC6A.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GPRC6A; G-protein coupled receptor family C group 6 member A; hGPRC6A; G-protein coupled receptor GPCR33; hGPCR33
Immunogen Species	Homo sapiens (Human)
Target Names	GPRC6A

Image



Western Blot analysis of Jurkat cells using GPRC6A Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.