




FGF13 Antibody

Product Code	CSB-PA008241
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92913
Immunogen	Synthesized peptide derived from the Internal region of Human FGF-13.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FGF13; FHF2; Fibroblast growth factor 13; FGF-13; Fibroblast growth factor homologous factor 2; FHF-2
Immunogen Species	Homo sapiens (Human)
Target Names	FGF13

Image



Western Blot analysis of JK cells using FGF-13 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.