




Phospho-FGFR1 (Y654) Antibody

Product Code	CSB-PA008243
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11362
Immunogen	Synthesized peptide derived from Human Flg around the phosphorylation site of Y654.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FGFR1; BFGFR; CEK; FGFBR; FLG; FLT2; HBGFR; Fibroblast growth factor receptor 1; FGFR-1; Basic fibroblast growth factor receptor 1; BFGFR; bFGF-R-1; Fms-like tyrosine kinase 2; FLT-2; N-sam; Proto-oncogene c-Fgr; CD antigen CD331
Immunogen Species	Homo sapiens (Human)
Target Names	FGFR1

Image



Western Blot analysis of HeLa cells using Phospho-Flg (Y654) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.