




Phospho-FGFR4 (Y642) Antibody

Product Code	CSB-PA008250
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22455
Immunogen	Synthesized peptide derived from Human FGFR-4 around the phosphorylation site of Y642.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FGFR4; JTK2; TKF; Fibroblast growth factor receptor 4; FGFR-4; CD antigen CD334
Immunogen Species	Homo sapiens (Human)
Target Names	FGFR4

Image



Western Blot analysis of JK cells using Phospho-FGFR-4 (Y642) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.