





RHOBTB3 Antibody

Product Code	CSB-PA008260
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94955
Immunogen	Synthesized peptide derived from the Internal region of Human RHOBTB3.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IHC:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RHOBTB3; KIAA0878; Rho-related BTB domain-containing protein 3
Immunogen Species	Homo sapiens (Human)
Target Names	RHOBTB3

Image



Western Blot analysis of K562 cells using RHOBTB3 Polyclonal Antibody



Western blot analysis of customer's lysis using RHOBTB3 antibody.