




ICK Antibody

Product Code	CSB-PA008268
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPZ9
Immunogen	Synthesized peptide derived from Human Intestinal Cell Kinase around the non-phosphorylation site of Y159.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ICK; KIAA0936; Serine/threonine-protein kinase ICK; Intestinal cell kinase; hICK; Laryngeal cancer kinase 2; LCK2; MAK-related kinase; MRK
Immunogen Species	Homo sapiens (Human)
Target Names	ICK

Image



Western Blot analysis of 3T3 cells using Intestinal Cell Kinase Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.