




RBM34 Antibody

Product Code	CSB-PA008308
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42696
Immunogen	Synthesized peptide derived from the C-terminal region of Human RBM34.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RBM34; KIAA0117; RNA-binding protein 34; RNA-binding motif protein 34
Immunogen Species	Homo sapiens (Human)
Target Names	RBM34

Image



Western Blot analysis of 3T3 cells using RBM34 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.