




MAST3 Antibody

Product Code	CSB-PA008310
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60307
Immunogen	Synthesized peptide derived from the C-terminal region of Human MAST3.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MAST3; KIAA0561; Microtubule-associated serine/threonine-protein kinase 3
Immunogen Species	Homo sapiens (Human)
Target Names	MAST3

Image



Western Blot analysis of RAT-HEART cells using MAST3 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.