




# FOX E1 Antibody

|                            |                                                                                                                                                                                                             |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA008315                                                                                                                                                                                                |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                               |
| <b>Uniprot No.</b>         | O00358                                                                                                                                                                                                      |
| <b>Immunogen</b>           | Synthesized peptide derived from the N-terminal region of Human FoxE1.                                                                                                                                      |
| <b>Raised In</b>           | Rabbit                                                                                                                                                                                                      |
| <b>Species Reactivity</b>  | Human,Mouse,Rat                                                                                                                                                                                             |
| <b>Tested Applications</b> | WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000                                                                                                                                             |
| <b>Form</b>                | Liquid                                                                                                                                                                                                      |
| <b>Conjugate</b>           | Non-conjugated                                                                                                                                                                                              |
| <b>Storage Buffer</b>      | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                     |
| <b>Purification Method</b> | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                       |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                         |
| <b>Alias</b>               | FOX E1; FKHL15; FOX E2; TITF2; TTF2; Forkhead box protein E1; Forkhead box protein E2; Forkhead-related protein FKHL15; HFKH4; HNF-3/fork head-like protein 5; HFKL5; Thyroid transcription factor 2; TTF-2 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                                                        |
| <b>Target Names</b>        | FOX E1                                                                                                                                                                                                      |

## Image



Western Blot analysis of HUVEC cells using FoxE1 Polyclonal Antibody

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.