



# FOXO1/FOXO3 Antibody

<b>Product Code</b>	CSB-PA008320
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q12778/O43524
<b>Immunogen</b>	Synthesized peptide derived from Human FoxO1/3 around the non-phosphorylation site of S322/S325.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Tested Applications</b>	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	FOXO1; FKHR; FOXO1A; Forkhead box protein O1; Forkhead box protein O1A; Forkhead in rhabdomyosarcoma; FOXO3; FKHL1; FOXO3A; Forkhead box protein O3; AF6q21 protein; Forkhead in rhabdomyosarcoma-like 1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	FOXO1/FOXO3
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.