



# Phospho-FOXO1/FOXO3/FOXO4 (T24/32) Antibody

<b>Product Code</b>	CSB-PA008321
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q12778/O43524/P98177
<b>Immunogen</b>	Synthesized peptide derived from Human FoxO1/3/4 around the phosphorylation site of T24/32.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	FOXO1; FKHR; FOXO1A; Forkhead box protein O1; Forkhead box protein O1A; Forkhead in rhabdomyosarcoma; FOXO3; FKHL1; FOXO3A; Forkhead box protein O3; AF6q21 protein; Forkhead in rhabdomyosarcoma-like 1; FOXO4; AFX; AFX1; MLLT7; Forkhead box
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	FOXO1/FOXO3/FOXO4
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.