




Phospho-FOXO3 (S253) Antibody

Product Code	CSB-PA008324
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43524
Immunogen	Synthesized peptide derived from Human FoxO3A around the phosphorylation site of S253.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FOXO3; FKHL1; FOXO3A; Forkhead box protein O3; AF6q21 protein; Forkhead in rhabdomyosarcoma-like 1
Immunogen Species	Homo sapiens (Human)
Target Names	FOXO3

Image



Western Blot analysis of HELA cells using Phospho-FoxO3A (S253) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.