





FLT1 Antibody

Product Code	CSB-PA008334
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17948
Immunogen	Synthesized peptide derived from Human Flt-1 around the non-phosphorylation site of Y1048.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	FLT1; FLT; FRT; VEGFR1; Vascular endothelial growth factor receptor 1; VEGFR-1; Fms-like tyrosine kinase 1; FLT-1; Tyrosine-protein kinase FRT; Tyrosine-protein kinase receptor FLT; FLT; Vascular permeability factor receptor
Immunogen Species	Homo sapiens (Human)
Target Names	FLT1

Image



Western blot analysis of KB SH-SY5Y 293T lysis using Flt-1 antibody.



Western blot analysis of HeLa KB SH-SY5Y 293T lysis using Flt-1 antibody.