



ZC3H4 Antibody

Product Code	CSB-PA008336
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPT8
Immunogen	Synthesized peptide derived from the Internal region of Human ZC3H4.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	IHC, ELISA; IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ZC3H4; C19orf7; KIAA1064; Zinc finger CCCH domain-containing protein 4
Immunogen Species	Homo sapiens (Human)
Target Names	ZC3H4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.